



VP8500 Video Cart

Like nothing else on the market, ConMed Linvatec introduces the next generation of video carts. VP8500 Video Cart provides the Operating Room with clean lines, new features, and added flexibility. Providing a platform for High Definition or Standard Definition equipment, the VP8500 is designed to be as forward looking as the components on its shelves.

- Monitor yolk for the VP4632 monitor
- Five adjustable 14" shelves
- HD-SDI outlet for auxiliary monitor
- Stainless Steel arms
- Built in isolation transformer and cord reel
- Lockable rear doors
- CO₂ Tank holder




CLEARLY TRUE HIGH DEFINITION™

True HD™



VP8500 Video Cart



CLEARLY TRUE HIGH DEFINITION™



Catalog Number

VP8500
VP8510
VP8520
VP8530
VP8540

Description

Video Cart
Foot Pedal Storage
Touch Screen Monitor Mount
Infusion Hanger
17" Wide Shelf


Specifications:

	Height (in/cm)	Width (in/cm)	Depth (in/cm)
Maximum:	79.3 / 201.4	22.5 / 57.2	22.5 / 57.2
Body:	55.3 / 140.4	22.5 / 57.2	22.5 / 57.2

Note: Dimensions listed are only for the cart; additional equipment or accessories may increase the overall dimensions.

Wheels:	Quantity 4, 5" twin casters (2 with locks)
Power outlets:	Quantity 10
Power Cord:	15 ft. with North American Hospital Grade Plug
Electrical Ratings	115 VAC, 60Hz, 10.6 Amps

Designed with today's busy ORs in mind, the VP8500 combines ease of use, a small footprint, and a high level of functionality



Built in isolation transformer and cord reel



Optional 17" wide shelf for shaver or pump consoles



HD-SDI output on the rear spine

CONMED CORPORATION PRODUCT AREAS:

ARTHROSCOPY • ELECTROSURGERY • ENDOSCOPY • ENDOSURGERY • GASTROENTEROLOGY • INTEGRATED SYSTEMS • PATIENT CARE • POWERED INSTRUMENTS • PULMONOLOGY



11311 Concept Boulevard
Largo, FL 33773-4908
(727) 392-6464
Customer Service: 1-800-237-0169
FAX: (727) 399-5256
International FAX: +1 (727) 397-4540
email: customer_service@linvatec.com
www.linvatec.com