

SHUTT[®] SERIES 210[™]

An addition to the time-honored Shutt Manual Instrument family




Endurance... Strength... Excellence...

GOING THE DISTANCE



SHUTT® SERIES 210™ MANUAL INSTRUMENTS



90° Rotary Forceps
 ■ Blunt Tip Right
 ■ Blunt Tip Left



Aggressor Forceps
 ■ 15° Up
 ■ 15° Left
 ■ 15° Right




Retrograde Forceps
 ■ 15° Right Bite



Suture Retrieval Forceps
 ■ Straight



Micro Scissors
 ■ Straight



Grasping Forceps
 ■ Straight with Ratchet

We listened to surgeon needs and hit the drawing boards to create this new and exciting line. The Shutt Series 210 family consists of nine manual instruments, all of which are 3.4mm in diameter and 210mm in working length, allowing for easy entry through cannulae and ample length for working in deep joints such as the hip.

Series 210 instruments feature the same barber-tang handle as their Shutt predecessors, offering a sensitive touch and a large thumb ring for long-term comfort. A special coating on the actuator guarantees a smooth cutting action while precision-matched, diamond-honed tips assure clean, close cuts. These new Shutt instruments feature a stronger tube and are laser welded at the shaft tip and base. The heavy duty shear pins stand up to the toughest tissue and the glass bead finish is easy to clean and gives each instrument a highly professional look and feel.

Shutt Series 210 – Truly going the distance.

ORDERING INFORMATION

Description	Cat. No.
Series 210, 90° Rotary Forceps, Right.....	S21ORFR
Series 210, 90° Rotary Forceps, Left	S21ORFL
Series 210, 15° Up, Aggressor Forceps	S210AFU
Series 210, 15° Aggressor Forceps, Left.....	S210AFL
Series 210, 15° Aggressor Forceps, Right.....	S210AFR
Series 210, 15° Up, Mantis Retrograde Forceps, Right Bite.....	S210RGF
Series 210, Suture Retrieval Forceps, Straight.....	S210SRF
Series 210, Micro Scissors, Straight	S210MSC
Series 210, Grasping Forceps, Straight with Ratchet.....	S210GSP

CONMED CORPORATION PRODUCT AREAS:

ARTHROSCOPY • ELECTROSURGERY • ENDOSCOPY • ENDOSURGERY • GASTROENTEROLOGY • INTEGRATED SYSTEMS • PATIENT CARE • POWERED INSTRUMENTS • PULMONOLOGY



11311 Concept Boulevard
 Largo, FL 33773-4908
 (727) 392-6464
 Customer Service: 1-800-237-0169
 FAX: (727) 399-5256
 International FAX: +1 (727) 397-4540
 email: customer_service@linvatec.com
 www.linvatec.com