





Reciprocating Action Shaver Blades

Superior design with the customer in mind!

- Reciprocates as it rotates
- Reduces clogging
- Promotes cleaner cutting without tearing
- Patent pending design produces the reciprocating action without additional hardware
- Now in 5 styles in 3.5 and 4.2mm diameters



ReAct[®] Reciprocating Action Shaver Blades



Teeth angles and geometry optimized for cutting, not ripping, of tissue and bone.

Shininess of outer tube is minimized to reduce glare and create a clearer field of view during surgery.



The Magnet Retainer creates unique cam cutting action.

Blade recognition technology allows rapid identification of blades and burs.

Hubs are color coded to denote different blade and bur styles.

Unique cam action design



Superior design with the customer in mind!





- Patent pending CAM design that produces the reciprocating action
- Patented manufacturing method for the XtraSharp[™] blades
- Sharp and aggressive
- Reduced occurrence of clogging
- Cleaner cut with less tissue tearing
- Offered in several window configurations

Ordering Information

ConMed Linvatec's patented XtraSharp[™] technology hones precision cutting edges which means superior cutting performance!

ConMed Linvatec's patent pending ReAct[™] (Reciprocating Action) technology allows for an efficient cutting motion that produces aggressive, low clogging resection of soft tissue.

For use with C9820, C9824, C9828, C9964, MC9820, MC9828, D9820, D9824, D9920 and D9924 Handpieces.

Description	Cat. No.
 Great White[®] The ultimate blade for aggressive soft tissue resection, including removal of the ACL stump. 3.5 mm Great White [®] 35GRW-RA-ZZ 4.2 mm Great White [®] 42GRW-RA-ZZ	
 Gator[®] Ideal for aggressive meniscal trimming, joint debridement, plica and synovium removal. 3.5 mm Gator [®] 35GTR-RA-ZZ 4.2 mm Gator [®] 42GTR-RA-ZZ	
 Full-Radius Resector Excellent for all general soft tissue debridement applications. 3.5 mm Full Radius Resector 35FRR-RA-ZZ 4.2 mm Full Radius Resector 42FRR-RA-ZZ	
Cuda[®] All purpose, aggressive blade, ideal for extensive tissue removal. 3.5 mm Cuda [®] 35CUD-RA-ZZ 4.2 mm Cuda [®] 42CUD-RA-ZZ	
 Tiger[®] Excellent for all general soft tissue debridement, especially meniscectomy and chondroplasty. 3.5 mm Tiger [®] 35TGR-RA-ZZ 4.2 mm Tiger [®] 42TGR-RA-ZZ	

CONMED CORPORATION PRODUCT AREAS:

ARTHROSCOPY • ELECTROSURGERY • ENDOSCOPY • ENDOSURGERY • GASTROENTEROLOGY • INTEGRATED SYSTEMS • PATIENT CARE • POWERED INSTRUMENTS • PULMONOLOGY



11311 Concept Boulevard
 Largo, FL 33773-4908
 (727) 392-6464
 Customer Service: 1-800-237-0169
 FAX: (727) 399-5256
 International FAX: +1 (727) 397-4540
 email: customer_service@linvatec.com
 www.linvatec.com