



# Dry-Doc®

Cannula

## SMOOTH & RIBBED CANNULA ALL IN ONE

- Robust hold in skin, muscle, and capsule
- Smooth surface for insertion and removal of cannula
- Hex cap interface provides secure interaction of cannula & obturator
- Double star seals prevent leaks
- Click lock ratcheting stopcock permits precise fluid control
- Semitranslucent color permits visual monitoring of instruments & sutures

# Dry-Doc<sup>®</sup> Cannula Leak Prevention System

## Description

	<b>Cat. No.</b>
8mm x 85mm Dry-Doc Cannula, translucent yellow .....	C7368
Packaged Sterile, 6/box with disposable obturators	
8mm x 75mm Dry-Doc Cannula, translucent red.....	C7367
Packaged Sterile, 6/box with disposable obturators	
7mm x 85mm Dry-Doc Cannula, translucent blue .....	C7360
Packaged Sterile, 6/box with disposable obturators	
5mm x 85mm Dry-Doc Cannula, translucent green .....	C7350
Packaged Sterile, 6/box with disposable obturators	

## Reusable Obturators

Reusable Obturator, anodized yellow for use with C7368.....	C7395
Reusable Obturator, anodized Red for use with C7367 .....	C7390
Reusable Obturator, anodized blue for use with C7360.....	C7385
Reusable Obturator, anodized green for use with C7350.....	C7380



Hex-cap interface provides secure interaction of cannula & obturator



Round distal tip prevents damage to joint capsule, suture, and Shuttle Relay



5mm x85mm

7mm x85mm

8mm x85mm

8mm x75mm



## CONMED CORPORATION PRODUCT AREAS:

ARTHROSCOPY • ELECTROSURGERY • ENDOSCOPY • ENDOSURGERY • GASTROENTEROLOGY • INTEGRATED SYSTEMS • PATIENT CARE • POWERED INSTRUMENTS • PULMONOLOGY



11311 Concept Boulevard  
Largo, FL 33773-4908  
(727) 392-6464

Customer Service: 1-800-237-0169

FAX: (727) 399-5256

International FAX: +1 (727) 397-4540

email: [customer\\_service@linvatec.com](mailto:customer_service@linvatec.com)

[www.linvatec.com](http://www.linvatec.com)