



Hex Flex™

Cannula

THE PLIABLE, BENDABLE, ULTRA DURABLE, CANNULA

- Two distal spiral threads provide a sturdy, yet atraumatic hold in capsule
- Smooth shaft for easy manipulation of cannula
- Compatible with ConMed Linvatec's color coordinated reusable cannulated obturators
- Click-lock ratcheting stopcock for meticulous fluid control
- Semitranslucent color for unobstructed visualization
- Double seals prevent fluid loss
- Hex cap interface provides easier insertion into joint capsule

Hex Flex™ Atraumatic "Soft Threaded" Cannula



Atraumatic Threads



Smooth Shaft



Hex Interface



5mm x 85mm

7mm x 85mm

8mm x 85mm



Description

Cat. No.

- | | |
|---|-------|
| 8mm x 85mm Hex Flex Cannula, translucent yellow
Packaged Sterile, 6/box with disposable obturators | C7480 |
| 7mm x 85mm Hex Flex Cannula, translucent blue
Packaged Sterile, 6/box with disposable obturators | C7465 |
| 5mm x 85mm Hex Flex Cannula, translucent green
Packaged Sterile, 6/box with disposable obturators | C7450 |

Reusable Obturators

- | | |
|---|-------|
| Reusable Obturator, anodized yellow for use with C7480..... | C7395 |
| Reusable Obturator, anodized blue for use with C7465..... | C7385 |
| Reusable Obturator, anodized green for use with C7450..... | C7380 |

CONMED CORPORATION PRODUCT AREAS:

ARTHROSCOPY • ELECTROSURGERY • ENDOSCOPY • ENDOSURGERY • GASTROENTEROLOGY • INTEGRATED SYSTEMS • PATIENT CARE • POWERED INSTRUMENTS • PULMONOLOGY



11311 Concept Boulevard
 Largo, FL 33773-4908
 (727) 392-6464
 Customer Service: 1-800-237-0169
 FAX: (727) 399-5256
 International FAX: +1 (727) 397-4540
 email: customer_service@linvatec.com
 www.linvatec.com